egger flex / AB

25/40/60 shore

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

Instructions for use

Highest wearing comfort due to soft-flexible fitting, multipurpose use. Ex-tremely tear resistant, long-term durable and easy processing. Light bodied, prolonged working time for mass-production in the lab, skin-friendly.

1. Mixing and dosing

Extrude material from the cartridge with the inject pistol (REF 25100 or 25120). Place cartridge into inject pistol and remove closure cap by rotation. Eject a small amount of material (bleeding) until equal flow is being extruded from both orifices. Insert mixing cannulas into guide grooves on cartridge and cannula and rotate clockwise to a stop. Dispenser is now ready to mix silicone in any amount needed. Store cartridge with mixing cannula attached for subsequent use.

2. Application in the lab – manuel

Prepare, wax or lacquer the impression for processing. All commercially available plasters and egger duplicating material (REF 40900 ff) can be used for the fabrication of the negative form (Fig. 1). Coat the plaster negative form with Isolat (REF 39900 ff). An insulation of gel forms is not necessary. The bubble-free injected material vulcanizes at room temperature without pressure application. If a handle is needed, we recommend using the Grip Stencil Set (REF 36899) according to the instructions for use. A plaster negative form with counter can be made for the fabrication of otoplastics, in order to shorten trimming: Fill one part of the brass flask with plaster and embed the insulated impression. After setting of the plaster deflask the impression. In the negative form of the impression pre-model the surface of the future earmold in wax. Insulate the plaster surface with a separating agent (e. g. dish washing liquid), mount the second part of the flask and fill it with plaster. After setting of the plaster open the flask. Scald both parts of the flask with boiling water and coat them with Isolat (REF 39900 ff). Extrude egger flex/AB from the cartridge (described under point 1) into the negative form and close flask with counter accurately (Fig. 2a). After vulcanization open the flask and remove the blank. In order to avoid airbubbles (e. g. with Aquapress REF 49800) curing can also be performed under pressure (Fig. 3).

3. Final trimming and varnishing

For molding and surface trimming use the special egger EF-cutters (for silicones) and trimmer (Fig. 4). If a handle is needed, we recommend using the Grip Stencil Set (REF 36899) according to the instructions for use. For smoothing the surface, use the air-drying egger silicone lacquers, Lacquer L (REF 37500 ff) and Lacquer L nano (with antibacterial effect, REF 37510 ff), Lacquer opaco (REF 37527) or the hot-vulcanizing Lacquer H/II (REF 39400 ff) according to their instructions for use (Fig. 5).

Important working hints

- Do not combine with condensation curing silicones.
- Cured silicone materials are chemically inert – spots on clothing should be avoided.
- Latex gloves and latex contaminated surfaces as well as cerumen, cremes and resins may inhibit the setting reaction of egger flex/AB (we recommend standard gloves made of nitrile rubber or polyethylene).
- The cartridges are filled in a vacuum, during quality control thoroughly tested and only released in air bubble free conditions.
- Since under certain environmental conditions (e. g. temperature, air pressure) air bubbles may be produced afterwards in the cartridge, curing of the material in a pressure-container is generally recommended.

Safety advice

The impression taking of the external auditory canal has to be executed by trained specialist staff only. The working instructions and precautionary measures have to be strictly observed. Non-observance might lead to irreparable damage of the ear or ear drum.

egger is not liable for any damage caused by improper application of the impression material.

Further information

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

Indications for use

Indirect fabrication of BTE and hearing protection earmolds, mass-production.

Technical Data

Mixed volume: 50 ml (cartridge)

Mixing ratio: 1:1

Product colors: available in various colors (see catalog)

Mixing time: omitted (Automix system)

Working time: approx. 3 min. 30 sec. (from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.)

Setting time: approx. 25–30 min. (from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.)

Final hardness: 25 / 40 / 60 shore A

Recovery from deformation: > 99.8%

Linear dimensional change: 0.2%

Application: at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity

Storage:

CAUTION:

Federal U.S. law restricts this device to sale by or on order of a licensed healthcare practitioner (or trained specialist personnel).

Please follow the product information and the material safety data sheet.