

Impression Syringe Standard



Technical Data Sheet

Intended use: Accessories for ear impressions

Composition:

Pos	chemical term
1	outer part: polycarbonate, transparent
2	inner part: polycarbonate-acrylonitrile butadiene styrene blend, opaque

Color Identification: outer part: clear transparent / pink transparent
 inner part: white opaque / pink opaque

Outer part:

Inner diameter (tip): approx. 2.5 mm

Outer diameter: approx. 16.5 mm

Length: approx. 130.0 mm

Filling volume: approx. 10 ml

Dreve Otoplastik GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 82909

otoplastik@dreve.de

www.dreve.com

These data result from measurements which are executed from time to time
in view of our quality management system.
This document is valid without any signature.

Rev. 3 / 04.2022
Page 1 of 1
TDS_Dreve Impression
Syringe Standard_GB