## Otoform® Xpand



## **Technical Data Sheet**

Product description: Addition vulcanizing ear impression silicone

Manufacturer: Dreve Otoplastik GmbH

Max-Planck-Straße 31 59423 Unna / Germany

Composition:

Pos	chemical term
1	vinyl containing polydimethylsiloxane
2	silicic acid
3	filling materials
4	pigments
5	platinum catalyst

Color: apricot

(corresponds to released sample in the house)

Processing time:  $100 \sec \pm 15 \sec$  at 23 °C

Setting time: 3 min  $\pm$  15 sec at 37 °C

Final hardness:  $20 \pm 2$  Shore A

Consistency: medium bodied type 2 (according to DIN EN ISO 4823)

Recovery from deformation: > 99 % (according to DIN EN ISO 4823)

Linear dimensional change: < 0.2 % (according to DIN EN ISO 4823)

Expansion of volume (after filling):  $3.5 \pm 1.5 \%$ 

Dreve Otoplastik GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 82909

otoplastik@dreve.de

www.dreve.com